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如遇停電，請將中間下方開關“ CPU & INSTRUMENT“ “ LAMP POWER & COOLING SYSTEM” 關閉

AVIV CD startup procedures

1. Open the valve on the liquid nitrogen dewar , check left pressure gauge indicate at 20 spi , the flow rate should be between 15-20 SCFH , waiting for 20 minutes to purge system.
2. startup the chiller at right side of the CD (ON/OFF button at chiller control panel)
3. Switch on the switch “ **LAMP POWER & COOLING SYSTEM** “ on POWER box , the LED “**CONTROL POWER ON**” will be on
4. wait until “**LAMP READY**” lit , press red button “ **IGNITE LAMP**” to ignite the lamp , “**LAMP ON**” lit indicate the light source was on .
5. Switch on the switch “ **CPU & INSTRUMENT**“ at POWER box , the LED “**CPU/INST PWR ON**” will on .
6. power on the PC at down left side
7. startup the CDS program
8. Display → Data browser → set “default data path “ at the bottom of the window -→ return

AVIV CD shutdown procedures

1. save your data and close CDS program
2. shutdown PC
3. Switch off the switch “ **CPU & INSTRUMENT** “ at **POWER** box
4. Switch off the switch “ **LAMP POWER & COOLING SYSTEM** “ on **POWER** box
5. shutdown the chiller at right side of the CD (**ON/OFF** button at chiller control panel)
6. waiting for 5 minutes and turn off(clockwise) the valve on liquid nitrogen tank

in case of power failure , please turn off both “ **CPU & INSTRUMENT** “ and “ **LAMP POWER & COOLING SYSTEM** “