# Malvernpanalytical Microcal PEAQ-ITC

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
Part	Description
1	Cell Cleaning Tool
2	Detergent bottle
3	Methanol bottle
4	Water bottle
5	Waste bottle
6	Overflow bottle

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
建議設定:  
**Reference power: 5 µcal/s**  
**Feedback: high**  
**Stir speed: 750 rpm**  
**Initial delay: 60 s**  
**Spacing time: 150 s**  
**Duration time: 4 s**  
**Number of injections: 19**  
**Injection volume: 2 µL**

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

Next

1. Choose Cleaning Method(s)
2. Insert Cell Cleaning Tool
3. Attach Fill Port Adapter
4. Move Pipette to Clean Location
5. Detach Fill Port Adapter
6. Remove Cell Cleaning Tool

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
**1 Load Cell**

Move the pipette out of the way (i.e. to the Clean Location).

Fill the loading syringe with 300  $\mu$ l of sample.

Slowly insert the loading syringe into the sample cell port, gently touch the cell bottom, and move up approximately 1 mm.

Slowly dispense approximately 150  $\mu$ l of sample.

Introduce a small volume, quickly, to dislodge any bubbles. Repeat this several times.

Slowly dispense the remaining sample while being careful not to introduce air bubbles.

Remove any excess sample in the cell port overflow cup using the loading syringe.

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**2 Attach Fill Port Adapter**

If the pipette is in the Clean Location, you must press the clamp's release lever.

Move the pipette to the Rest Location.

Align the hole in the pipette's housing to the hole in the pipette's rotating assembly.

Insert the fill port adapter. A soft click should be felt.

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**3 Move Pipette to Load Location**

Load approximately 60  $\mu$ l of titrant in one of the supplied microcentrifuge tubes.  
 Ensure the microcentrifuge tube has its lid properly situated in the keyed Load Location.  
 Move the pipette to the Load Location.  
 Click Next.

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**5 Move Pipette to Cell**

If the cell is loaded, move the pipette into the cell.  
 Be sure the pipette is firmly seated in the cell port.  
 Now you may start your experiment.  
 Click Done

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**Run**

Experiment Information


[Syr] (M) 5.00e-3  
 [Cell] (M) 400e-6  
 Comment #cell-EDTA, #syringe-CaCl2

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