

**Transgenic Core Facility**  
**Institute of Molecular Biology, Academia Sinica**  
**2789-9312, 2652-1438**

Case #: \_\_\_\_\_

**Gene Targeting Construct DNA Submission Form**

<b>Construct Name</b>	
<b>DNA Preparation Method</b>	
<b>DNA Preparation Date</b>	
<b>Linearized DNA Quantity</b>	
<b>Linearization Date</b>	
<b>DNA Submission Date</b>	
<b>I submitted _____<math>\mu</math>g DNA in _____<math>\mu</math>l of digestion reaction sol'n.</b>	
<b>The gel photo I pasted below is taken on _____(date)</b>	
<p><b>Note: Please make sure you satisfy the following requirements</b></p> <p><input type="checkbox"/> Gel electrophoresis should be clear and all the fragments should be fully separated</p> <p><input type="checkbox"/> Gel photo should be large and the following info should be included:</p> <ul style="list-style-type: none"> <li>a. MW marker</li> <li>b. Uncut circular plasmid</li> <li>c. Linearized fragment (100ng)</li> <li>d. Linearized fragment (500ng)</li> </ul> <p><input type="checkbox"/> All fragments should be clearly indicated by size and name</p>	<p style="text-align: center;"><b>(Please paste your gel photo here)</b></p>
<b>Requester: _____, Date: _____</b>	